

DIP MOLDING

FULL-SERVICE CAPABILITIES

At Rose Medical, we specialize in dip molding of custom, compliant balloons. Whether you're needing rapid prototypes or large-volume production, our team of experts ensure precision, quality, and consistency at every stage. By meticulously controlling each step of the process in-house, we are able to meet tight tolerances and the unique requirements of your device, on time and to specifications.

CUSTOM CAPABILITIES

Our most common materials include silicone and urethane, with typical wall thicknesses from 0.002" to 0.020".

Expertly manufactured with smooth or textured surfaces for enhanced performance. Along with our overmolding capabilities, we are ready to meet and exceed your devices unique requirements.



ISO 13485 Certified



FDA Registered



Class 7 Cleanrooms